

## Taq DNA Polymerase (with dNTPs), Economy

02-001 200 U (5 U/ µl), 02-001-5 5 x 200 U (5 U/ µl)

**Storage:** Ship at  $4^{\circ}$ C or  $-20^{\circ}$ C and store at  $-20^{\circ}$ C.

## Concentration: 5 units/ µl

\* Note: One unit is defined as the amount of enzyme that can incorporate 10 nmols of total dNTPs into an acid-insoluble material in 30 minutes at 74°C when activated salmon sperm DNA is used as template/primer.

**Storage Buffer:** 20mM Tris-HCl (pH 8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.5% Tween 20, 0.5% Igepal CA-630

**Supplied Reagents:** 10 x Standard Buffer (*Taq*): 100mM Tris-HCl (pH 8.3), 500mM KCl, 15mM MgCl<sub>2</sub>, 2.5mM (each) dNTPs

## Applications:

- 1) High-throughput PCR
- 2) Colony PCR
- 3) Incorporation of dUTP, dITP, and fluorescence-labeled nucleotides
- 4) Primer extension
- 5) Addition of a single nucleotide (adenosine) at the 3'-blunt ends

**Background:** *Thermus aquaticus* DNA polymerase (*Taq* DNA polymerase) was expressed in *E. coli* in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and the MW is 94 kDa. This enzyme is suitable for PCR reactions; capable of amplifying DNA with various primers.

Quality Assurance: Greater than 95% purity as determined by SDS-PAGE (CBB staining) (Fig.1)

The absence of endonucleases and exonucleases was confirmed.

**PCR Test:** Good amplification result was obtained in PCR reaction using  $\lambda$ DNA as a template (Fig.2).

## Related product:

# <u>02-021</u> Pfu DNA polymerase (+dNTPs), Economy

# <u>02-031</u> Pfu DNA polymerase (-dNTPs), Economy

<u>General composition of PCR reaction mixture (total 50 µl)</u>		
Taq DNA polymerase (5 units/ul)		*0.25 µl
10 x Standard Buffer ( <i>Taq</i> )		5 µl
2.5mM (each) dNTPs		4 µl
Template		<500 ng
Primer 1	$0.2{\sim}1.0\mathrm{p}$	ıM (final conc.)
Primer 2	$0.2{\sim}1.0\mathrm{p}$	ıM (final conc.)
Sterile distilled water		up to 50 µl
*Use of excess amount is not recommended.		



25 cycles





Fig.1 SDS-PAGE of TaqDNA polymerase



 $\begin{array}{ll} M:marker,\\ lane 1:2 \ kb, & lane 2:4 \ kb,\\ lane 3:6 \ kb, & lane 4:8 \ kb. \end{array}$ 

