

## cDNA Library, *S. cerevisiae*, Log Phase

02-701 500 ng

This cDNA library (plasmid DNA) is constructed from *Saccharomyces cerevisiae*, strain S288C-derived poly(A)<sup>+</sup> RNA at the log phase by the Linker-Primer method (Ref.1) by Prof. H. Nojima of Osaka University. This library is unidirectionally cloned by using the oligo (dT)<sub>18</sub> linker primer which contains the restriction enzyme site of *Not* I, and *Bam*HI (*Bgl* II)-*Sma* I adaptor.

The pLZ3 vector (shown below) used in this library can not replicate in *S. cerevisiae* but contains pUCori for replication in *E. coli*

### Application

PCR screening of known or unknown gene: Prepare the primers for the known or unknown gene (cDNA) and amplify the gene by PCR from this library followed by cloning to an appropriate vector.

Standard amplifying conditions: 35 cycles of PCR reactions using 10-100 ng of cDNA as a template. (Change the quantity of template and the number of cycles depending on the expression rate of mRNA of the objective gene.)

### Specification

**Quantity:** 500 ng (40 ng/ul, 13ul) in 10 mM Tris-HCl-1mM EDTA (pH 7.5)

**Quality:** 1) Number of independent clones:  $3.6 \times 10^6$   
2) Average insert size : longer than 1 kb

**Storage:** -20°C

**References:** Construction of this library is described in Supplementary data of Ref.3


1. Kobori M *et al* "Large scale isolation of osteoclast-specific genes by an improved method involving the preparation of a subtracted cDNA library." *Genes Cells* **3**: 459-475 (1998) PMID: [9753427](#)
2. Tanaka S and Nojima H "Nik1: a Nim1-like protein kinase of *S. cerevisiae* interacts with the Cdc28 complex and regulates cell cycle progression." *Genes Cells* **1**, 905-921 (1996) PMID: [9077450](#)
3. Tougan T, Okuzaki D, Nojima H. Chum-RNA allows preparation of a high-quality cDNA library from a single-cell quantity of mRNA without PCR amplification. *Nucleic. Acids Res.*, 36(15):e92, (2008) PMID:[18603591](#)

### Note

- \* This library is to be used only by the purchaser. It is not allowed to amplify and transfer the library to a third person.
- \* Related products: human tissue specific cDNA libraries and cDNA libraries of model organisms (See [HP](#)).

to be continued...

Fig. Structure of pLZ3 and the restriction sites.



; MCS(pLZ3)

```

CpoI(3)  SauI(b) MluI(5)                AatII(3) BglII(5) AscI(5)  BalI(b)
PstI(3) SacI(3)  ApaI(3)  -----
SseI(3) -----                T7 Promoter  EcoRI(5)  XbaI(5)  AflIII(5)  BstXI(5)
-----
NNNCTGCA  CCTGCAGGAGCTCGGACCGGGCCCTTAGGACGGTAATACGACTCAGTATAGGGAATTCGACGCTAGATCTTAAGGGCGGCCAAGGGGTTGGCCA
NNNG  ACGTGGACGCTCCTCGAGCCTGGCCCGGGAATCTGCGCATTATGCTGAGTGATATCCCTTAAGCTGCAGATCTAGAATTCGCGCGGGTCCCAACCGGT

BstEII(5)
-----
SnaBI(b)  DraIII(3)  NheI(5)  -----  SceI(3)  NotI(5)  T3 promoter  SwaI(3)  SphI(5)  NruI(b)  PacI(3)  SacII(3)  SacI(3)
-----
CGTGGTAACCACGGGGTGGCTAGCTAGGGATAACAGGGTAATATAGCGGCCGCCCTTAGTGAGGGTTAATTTAAATCGTACGTCGGGATTAATTAACCGCGGTGGAGCT  CAAT
GCACCATTTGGTCCCGCCACGATCGATCCCTATTGTCCCATATATCGCCGGCGGAAATCACTCCCAATTAATTTAGCATGCAGCGCTAATTAATTTGGCGCCACC  TCGACTTA

TCGCCCTATAGTGAGTCGTATTA -3'
AGCGGGATATCACTCAGCATAAT -5'

```