

DNA polymerase kappa, functional

10-105 50 ug

Mammalian **DNA polymerase** κ , a member of the UmuC/DinB nucleotidyl transferase superfamily, has been implicated in spontaneous mutagenesis (1). Human **DNA polymerase** κ copies undamaged DNA with average single-base subsitution and deletion error rates of 7 x 10⁻³ and 2 x 10⁻³, respectively. These error rates are high when compared to those of most other DNA polymerases (2). **DNA polymerase** κ has important role in the mutagenic bypass of certain types of DNA lesions (3).

This product was over-expressed as a recombinant protein in *E. coli* with a plasmid carrying a *C*⁴terminal histidine-tagged human **DNA polymerase** κ (1-560 aa), and highly purified by several steps of chromatography (2). The product is catalytically active and its molecular weight is 65 kD (Fig 1). Activity of this product has been confirmed by a user researcher even if it was diluted 8,000-fold.

Application

1) Functional studies on translesion DNA synthesis by DNA polymerase kappa 2) Antigen for western blotting and ELISA by anti-Pol kappa antibody Purity: Over 90% by SDS-PAGE (CBB staining) Protein concentration: 3.2 mg/ml as measured by BCA method Form: 0.2 M NaCl, 10 mM sodium phosphate buffer (pH 7.0), 50% glycerol Storage: Sent at 4°C or at -20°C. Upon arrival, centrifuge briefly and store at -20° C or at -80° C for longer storage. Data Link: UniProt KB/Swiss-Prot __Q9UBT6_ (POLK_HUMAN) **References**: This product has been used in the following publications 1. Ohashi E et al (2000) Fidelity and processivity of DNA synthesis by DNA polymerase kappa, the product of the human DINB1 gene. J Biol Chem 275: 39678-39684 (2000) PMID: <u>11006276</u> 2. Ohashi E et al (2000) Error-prone bypass of certain DNA lesions Fig.1 by the human DNA polymerase kappa" Genes Dev 14: 1589-1594 (2000) PMID: _10887153_ 3. Jałoszyński P. et al (2005) Error prone and inefficient replication across 8-hydroxyguanine (8-oxoguanine) in human and

.mouse _*ras_* gene. _Genes Cells. 10:543-50. PMID: _<u>15938713</u>_



Fig. 1.DNA polymerizing activity of Pol kappa (1-560 aa): 50 mM Tris-HCl (pH 8.0), 2 mM MgCl₂, 1 mM DTT, 0.1 mg/ml BSA, 0.1 mM dNTP, 100 nM primer/template (13-mer/30-mer) DNA polymerase κ 13 - 50 ng/25 ul assay, at 37°C, 10 min (above figure)

Fig.2 polyacrylamide gel electrophoresis of DNA polymerase κ protein

50

BioAcademia,Inc. Tel. 81-6-6877-2335 Fax. 81-6-6877-2336 info@bioacademia.co.jp http://www.bioacademia.co.jp/en/