Whatman™ Mini-UniPrep G2 Syringeless Filters

Product Information sheet

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.

Do not use internally or externally in humans or animals.

Introduction

Mini-UniPrep™ G2 Syringeless Filters are preassembled filtration units for removal of particulates from low volume liquid samples.

They are primarily designed for sample preparation prior to High Performance Liquid Chromatography (HPLC)/Ultra High Performance Liquid Chromatography (UHPLC). Mini-UniPrep G2 replaces syringe filters, syringes and autosampler vials/septa/caps with a single, disposable unit.


Each Mini-UniPrep G2 unit consists of two parts: a glass chamber and a plunger with an inner glass reservoir. The design incorporates a filtration membrane on one end of the plunger and a pre-attached septum cap on the other. By pressing the plunger through a liquid sample contained in the chamber, positive pressure forces the filtrate into the inner reservoir of the plunger. Air escapes through the vent hole until the evaporation seal is engaged, which then provides an airtight seal. The assembled unit can be used with any auto-sampler that takes standard 12 mm x 32 mm profile vials.


The Mini-UniPrep G2 is constructed with borosilicate glass, designed to maintain the purity of your sample by reducing the exposure of sample solvent to plastic components, which can potentially contaminate the sample.

- Before filtration, the sample is contained within a borosilicate glass outer chamber only.

- During filtration, the filtered sample is directed into the inner surface of the plunger and is then deflected into an inner reservoir also comprised of borosilicate glass.
- After filtration, the filtrate aliquot available for downstream analysis is contained entirely within the inner glass insert.

Therefore, during normal use, the sample will contact plastic components only during the very short time covering the actual act of filtration (i.e. only during compression of the plunger).

Safety Precaution

ALWAYS USE A MINI-UNIPREP G2 COMPRESSOR DEVICE.

Do not attempt to compress the device by hand; there is a risk of breakage of the glass chamber which could result in personal injury if a compressor is not used. Mini-UniPrep starter packs contain 100 Mini-UniPrep G2 syringeless filters and a Hand Compressor. If a Hand Compressor is not available do not use the Mini-UniPrep G2 Syringeless filters until one has been obtained*.

Do not use any other device to compress the Mini-UniPrep G2 filters including the six-position compressor (catalogue no. CR0000006) designed to be used with the original Mini-UniPrep filters.

Do not overfill the Mini-UniPrep G2 or compress the device at an angle before the airtight seal is made. Liquid could escape from the air vent holes, resulting in the loss of sample or personal injury.

* A Hand Compressor (catalogue no. MUPG2PWC1) is available to order separately.

See www.whatman.com for further details.


Instructions for Use

Before you start, remember to select the Mini-UniPrep G2 Syringeless Filter that is compatible with your sample (see Chemical Compatibility Guide in appendix 1).

Step 1 Insert a glass chamber into the Hand Compressor base.	
Step 2 Dispense the sample into the glass chamber taking care not to overfill. The Mini-UniPrep G2 glass chamber has a maximum capacity of 500 µl indicated by a printed "Fill Line" on the glass chamber. Minimum sample volume that may be added to the glass chamber is 220 µl in order to collect 50 µl in the glass insert.	
Step 3 Place a plunger into the neck of the glass chamber. Do not attempt to depress the plunger further than the neck of the glass chamber.	
Step 4 Place the Hand Compressor lid over the device such that the recess on its inner face fits over the septum cap of the plunger and the alignment guide is orientated with the slot in the base.	
Step 5 Steady the Hand Compressor base with one hand on a flat surface whilst using the other hand to press down on the Hand Compressor lid until the plunger is fully depressed into the chamber.	

Step 6

Lift the Hand Compressor lid from the base and remove the Mini-UniPrep G2.

The Mini-UniPrep G2 prepared sample is now ready for analysis. Process in either of two ways:

- A: Place the Mini-UniPrep G2 into an autosampler that takes standard 12 mm x 32 mm profile vials for automated sample injection into your instrument. Before your run, ensure that the autosampler needle height is set to a minimum height of 4 mm, to avoid contact with the bottom of the inner glass insert.
- B: Pierce the septum and draw the filtered sample into a syringe for manual injection into your instrument.

For optimum performance, the assembled and filled Mini-UniPrep G2 device should be kept in an upright position and not be subjected to excessive agitation to prevent the sample being spilled out of the protective inner glass storage reservoir.

Ordering Information

Mini-UniPrep G2 Syringeless Filters are available in different variants such as different membrane types, amber colored version for light sensitive samples and slit septum caps.

Please visit www.whatman.com for more information.

Appendix

Chemical Compatibility of Membrane

Solvent	GMF	NYL	PES	PP	PTFE	PVDF	RC	DpPP
Acetic Acid 5% +	R	R	R	R	R	R	R	R
Acetic Acid, Glacial	R	LR	R	R	R	R	NR	R
Acetone	R	R	NR	R	R	NR	R	R
Acetonitrile	R	R	NR	R	R	R	R	R
Ammonia, 6N	LR	R	R	R	R	LR	R	R
Amyl Acetate	R	R	LR	R	R	LR	R	R
Amyl Alcohol	R	R	NR	R	R	R	R	R
Benzene *	R	LR	R	LR	R	R	R	LR
Benzyl Alcohol *	R	LR	NR	R	R	R	R	R
Boric Acid	R	LR	+	R	R	R	R	R
Butyl Alcohol	R	R	R	R	R	R	R	R

Solvent	GMF	NYL	PES	PP	PTFE	PVDF	RC	DpPP
Butyl Chloride *	R	NR	+	NR	R	R	R	+
Carbon								
Tetrachloride *	R	LR	NR	LR	R	R	R	LR
Chloroform *	R	NR	NR	LR	R	R	R	LR
Chlorobenzene	R	NR	NR	LR	R	R	R	+
Citric Acid	R	LR	R	R	R	R	R	+
Cresol	R	NR	NR	R	R	NR	R	R
Cyclohexane	R	R	R	R	R	R	R	R
Cyclohexanone	R	NR	NR	R	R	R	R	R
Diethyl Acetamide	R	R	+	R	R	NR	R	R
Dimethyl								
Formamide	R	R	NR	R	R	NR	LR	R
Dioxane	R	R	LR	R	R	LR	R	R
DMSO	R	R	NR	R	R	LR	LR	R
Ethanol	R	R	R	R	R	R	R	R
Ethers	R	R	R	R	R	LR	R	R
Ethyl Acetate	R	R	NR	R	R	LR	R	R
Ethylene Glycol	R	R	R	R	R	R	R	R
Formaldehyde	R	R	R	R	R	R	R	R
Formic Acid	R	NR	R	R	R	R	LR	R
Freon TF	R	R	R	R	R	R	+	R
Hexane	R	R	R	R	R	R	R	R
Hydrochloric								
Acid (Conc)	R	NR	R	LR	R	R	NR	LR
Hydrofluoric Acid	NR	NR	+	LR	R	R	NR	LR
Isobutyl Alcohol	R	R	+	R	R	R	R	R
Isopropyl Acetate	R	R	+	R	R	R	R	R
Methanol	R	R	R	R	R	R	R	R
Methyl Ethyl								
Ketone	R	R	NR	R	R	NR	R	R
Methylene								
Chloride *	R	NR	NR	LR	R	R	R	LR
Nitric Acid (Conc)	R	NR	NR	NR	R	R	NR	NR
Nitric Acid, 6N	R	NR	LR	LR	R	R	LR	LR
Nitrobenzene *	R	LR	NR	R	R	R	R	R
Pentane	R	R	R	R	R	R	R	LR
Perchloro Ethylene	R	R	NR	R	R	R	R	LR
Phenol (0.5%)	R	R	NR	R	R	R	R	R
Pyridine	R	LR	NR	R	R	R	R	R
Sodium								
Hydroxide, 6N	NR	LR	R	R	R	R	NR	R
Sulfuric Acid (Conc)	R	NR	NR	NR	R	R	NR	R
Tetrahydrofuran	R	R	NR	LR	R	R	R	LR
Toluene *	R	LR	NR	LR	R	R	R	LR
Trichloroethane *	R	LR	NR	R	R	R	R	R
Trichloroethylene *	R	NR	NR	R	R	R	R	R
Water	R	R	R	R	R	R	R	R
Xylene *	R	LR	LR	LR	R	R	R	LR

Note: R = Resistant; LR = Limited Resistance; NR = Non Resistant; + = Insufficient Data; * = Short term resistance of housing.

Material Abbreviations

GMF = Glass Microfibre

NYL = Nylon

PES = Polyethersulfone

PP = Polypropylene

PTFE = Polytetrafluoroethylene

PVDF = Polyvinylidene Difluoride

RC = Regenerated Cellulose

DpPP = Polypropylene Depth Filter

Membrane Guide

Type	Applications
Glass Microfiber (GMF)	High loading capacity, binder-free, glass microfiber depth filter medium.
Nylon (NYL)	Membrane for aqueous and organic samples within a pH range of 3-10.
Polyethersulfone (PES)	For aqueous-based samples. Good flow rates.
Polypropylene (Dp PP)	Aqueous and organic samples; high loading capacity depth filter.
Polypropylene (PP)	For solvent-based samples, low water breakthrough values.
Polytetrafluoro-ethylene (PTFE)	For aggressive samples. Great chemical compatibility.
Polyvinylidene difluoride (PVDF) (Hydrophilic)	Low nonspecific protein binding membrane for samples in aqueous solutions or organic solvents.
Regenerated Cellulose (RC)	Suitable for use with either aqueous or organic solvents. Very low non-specific protein binding membrane.

Legal.

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

Whatman and Mini-UniPrep are trademarks of GE Healthcare companies.

Patents pending

© 2012 General Electric Company – All rights reserved.

First published May 2012.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit www.gelifesciences.com/contact

GE Healthcare UK Limited

Amersham Place
Little Chalfont, Buckinghamshire,
HP7 9NA, UK

<http://www.gelifesciences.com>

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala,
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg,
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo 169-0073,
Japan

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© 2012 General Electric Company - All rights reserved.