无菌采样阀 Aseptic Sampling Valve

Working Principle 工作原理

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- The Aseptic Sampling Valve is designed for sterilization before and after each sampling.
- Open valve: Sampling--The manual valve is opened by rotating the handle (adjustable valve) or by activating the lever(on/off valve).The stem and the membrane are then retracted, allowing liquid to pass.
- Closed valve: Sterilizing--The built-in spring closes the valve and keeps the channel open between the hoses though reverse process in order to sterilize. If with steam, a small pressure relief valve (optional) on the outlet is recommended. If using steam a small pressure relief valve (optional) on the outlet is recommended.
- Standard Design:The Aseptic Sampling valve consists of three parts, a valve body, a handle and a membrane. The rubber membrane is placed on the stem of the valve and works as a stretchable plug.
- Connection:Thread,Tank welding,Pipe welding,Clamp
- Handle:(1) Turning handle (2) Removable key
- Application: The Aseptic Sampling Valve is applied to beer,milk,beverage,fruit juice,pharmacy etc industries.

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



Ordering: It is recommended that the membrane be changed every 200 sterilizations.

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Working Principle Diagram

开启阀门, 取样
Open valve, Sampling.

关闭阀门, 消毒
Closed valve, Sterilizing.

序号 No.	名称 Description	材质 Material
1	阀体 Valve body	304 316L
2	密闭套 Closed pad	EPDM
3	阀杆 Valve shaft	304 316L
4	出口外套 Outlet outer sleeve	硅橡胶 Silicone rubber
5	弹簧 Spring	304
6	弹簧套 Spring sleeve	304 316L
7	螺母 Nut	304
8	手柄 Handle	304
9	手柄定位套 Handle locating sleeve	304
10	垫片 Pad	304
11	内六角螺钉 Inner socket blot	304
12	弹垫 Spring pad	65Mn
13	垫圈 Gasket	PTFE