
# **AppliedPhotophysics** More time for science

## Introduction

Understanding or preventing degradation is essential when developing a complex biotherapeutic. This study uses a novel approach to HOS analysis — quantitative circular dichroism — to detect and quantify statistically significant changes in HOS of IgG1 for important degradation pathways.

Use of a Chirascan-auto qCD system provided a rapid, highly productive, auditable method for HOS comparisons with:

- High turnaround of multiple samples
- Unattended, reproducible operation during spectral acquisition
- Objective quantification of subtle differences not discernible by eye
- Elimination of subjective visual comparisons and human error



4°C for seven of	days, se	erved as a	a referen


Aliquot 6

# Quantitative Circular Dichroism (<sub>q</sub>CD) for the Comparison of the Higher-Order Structures of Monoclonal Antibodies

Doug Marshall (doug.marshall@photophysics.com)<sup>1</sup>, <u>James Law</u><sup>1</sup>, Tom Hampson<sup>1</sup>, Dariusz Sliwa<sup>2</sup>, Lindsay Cole<sup>1</sup>, Li Luo<sup>2</sup> and David Gregson<sup>1</sup> <sup>1</sup>Applied Photophysics Limited 21 Mole Business Park, Leatherhead, Surrey, United Kingdom, KT22 7BA <sup>2</sup>Applied Photophysics Inc, 100 Cummings Center, Suite 440C, Beverly, MA 01915

### Experimental set-up prior to unattended operation

- Five aliquots from each treatment (#1-5) were presented on a 96-well plate.
- Aliquots of buffer were included for each sample (B).
- Three repeat Near UV CD and absorbance spectra were measured for each well using a 5 mm pathlength quartz flow cell. An equivalent Far UV experiment was performed at 0.1 mm pathlength.
- 180 CD and 180 absorbance spectra were measured in unattended operation for 18 hours (Near UV) or for 12




## Equipment

Samples were run and analyzed on a Chirascan-auto qCD. This fully automated system enabled repeat measurements of multiple samples thereby increasing signal: noise and enabling the generation of quantitative, statistically significant data.

Proprietary DichOS<sup>™</sup> calibration, standardized operating procedures and comprehensive statistical analysis combine to provide a novel, information-rich approach to biophysical characterization of proteins and peptides by circular dichroism.



Chirascan-auto qCD is the highest sample turnaround system in the Chirascan<sup>™</sup> platform







### Objective statistical analysis provided quantitative results

Two well-established statistical methods<sup>1, 2</sup>, correlation coefficient (Pearson coefficient) and weighted spectrum difference (WSD), were applied:

WS

$$r = \frac{\sum_{i=1}^{n} (y_i - \overline{y})(x_i - \overline{x})}{\sqrt{\sum_{i=1}^{n} (x_i - \overline{x})^2 - \sum_{i=1}^{n} (y_i - \overline{y})^2}}$$

$$SD = \sqrt{\sum_{i=1}^{n} \left[ \frac{1}{n} \left( \frac{|x_i|}{|x_i|_{ave}} \right) (x_i - y_i)^2 \right]}$$

Green = Similar ( $2\sigma$  significance)

Red = Dissimilar ( $2\sigma$  significance)

x<sub>i</sub> = signal from reference spectrum

 $y_i$  = signal from sample spectrum

n = number of data points

Characteristics: Pearson: 0 < r < 1, 1 = identical WSD > 0, 0 = identical

Similarity scores were calculated from a reference set of spectra and from sets of sample spectra. A Student's t-test was used to calculate the probability that the two populations of similarity scores were representative of the same population.

If p < 0.05, then difference is significant at the 95% confidence level.

#### Analysis of Near UV CD Da-

	Condition	Scores for comparison with control (aliquot 6)			
Sample		Correlation coefficient		Weighted spectral difference	
		Similarity	p-value (t-test)	Similarity	p-value (t-test)
Aliquot 1	pH8.5 and at 40°C for 1 week	0.99995	0.024	0.074	0.001
Aliquot 2	50 mM glucose at 20°C for 24 hours	0.99998	0.312	0.032	0.524
Aliquot 3	2 M glucose at 40°C for 1 week	0.99993	<0.001	0.050	0.011
Aliquot 4	0.3% H <sub>2</sub> O <sub>2</sub> at 20°C for 3 hours	0.99994	<0.001	0.180	<0.001
Aliquot 5	Control	0.99998	0.923	0.027	0.807

#### Analysis of Far UV CD Data:

Sample	Condition	Scores for comparison with control (aliquot 6)					
		Correlation coefficient		Weighted spectral difference			
		Similarity	p-value (t-test)	Similarity	p-value (t-test)		
Aliquot 1	pH8.5 and at 40°C for 1 week	0.99983	0.289	0.090	0.236		
Aliquot 2	50 mM glucose at 20°C for 24 hours	0.99993	0.103	0.052	0.193		
Aliquot 3	2 M glucose at 40°C for 1 week	0.99987	0.28	0.079	0.236		
Aliquot 4	0.3% H <sub>2</sub> O <sub>2</sub> at 20°C for 3 hours	0.99976	0.638	0.134	0.103		
Aliquot 5	Control	0.99987	0.418	0.076	0.358		

## Conclusions

- The affects of all but one of the treatments designed to induce deamidation/ isomerization, glycation and oxidation in IgG1 were detected in the Near UV CD datasets and shown to be significant at the  $2\sigma$  confidence interval.
- Significant differences in the equivalent Far UV datasets were not detected.
- Together, these results suggest that treatments caused changes to the local environment of aromatic side chains (tertiary structure), but did not cause changes to the secondary structure content.
- The ability to generate truly quantitative data will substantially strengthen the role of CD analysis throughout biotherapeutic development programs and help fulfil regulatory demands.

1. Dinh, N.N. et al. Quantitative spectral comparison by weighted spectral difference for protein higher order structure confirmation, Analytical Biochemistry 464 (2014), 60-62. 2. Teska, B.M. et al. Comparison of quantitative spectral similarity analysis methods for protein higher-order structure confirmation, Analytical Biochemistry, 434.1 (2013), 153-165