

Data Sheet

FluoTag-X2 anti-Guinea Pig Fluorophore-conjugated

Fluorophore	Catalog No. (20µl)	Catalog No. (200µl)
Atto488	N0602-At488-S	N0602-At488-L
Atto542	N0602-At542-S	N0602-At542-L
Atto565	N0602-At565-S	N0602-At565-L
Atto647N	N0602-At647N-S	N0602-At647N-L
Abberior-Star 580	N0602-Ab580-S	N0602-Ab580-L
Abberior-Star 635P	N0602-Ab635P-S	N0602-Ab635P-L
Sulfo-Cyanine 3	N0602-SC3-S	N0602-SC3-L
Sulfo-Cyanine 5	N0602-SC5-S	N0602-SC5-L

Identity Camelid sdAb anti-Guinea Pig

Content 5 µM fluorescently labeled sdAb in buffered saline,
50% glycerol, 0.09% sodium azide.

Storage Up to 3 months: -20 °C
Up to 12 months: -80 °C or below
Protect from light!

Applications IF / ICC Recommended dilution 1:500

Specificity Recognizes specifically Guinea Pig IgGs
Does not cross-react with antibodies from species
like: mouse, rabbit, rat, chicken or goat.

For more information please visit our web page at www.nano-tag.com

To be used *in vitro* / for research only.
Non-toxic, non-hazardous, non-infectious.