

Anti-HIV-1 Protease Antibody, rabbit serum

65-018 100 μl,

Background: HIV-1 protease is the aspartyl protease that mediates proteolytic cleavages of Gag and Gag-Pol polyproteins during or shortly after the release of the virion from the plasma membrane. Cleavages take place as an ordered, step-wise cascade to yield mature proteins. This process is called maturation. Displays maximal activity during the budding process just prior to particle release from the cell. Also cleaves Nef and Vif, probably concomitantly with viral structural proteins on maturation of virus particles. Hydrolyzes host EIF4GI and PABP1 in order to shut off the capped cellular mRNA translation. The resulting inhibition of cellular protein synthesis serves to ensure maximal viral gene expression and to evade host immune response.

Applications

- 1. Western blotting (1/2,000)
- 2. Dot blotting (1/2,000)
- 3. Immunoprecipitation (1/200)
- 4. Inhibition of HIV-1 Protease activity (assay dependent)
- 5. ELISA. (1/10,000)
 - Other applications have not been tested.

Immunogen: Full-size functional recombinant HIV-1 protease expressed and purified from E. coli

Reactivity: Reacts with HIV-1 protease of all substrains

Form: Whole antiserum added with 0.09% sodium azide

Shipping and Storage: Sent at 4° C or at -20° C and stored at -20° C. When necessary, freezing (with liquid or ethanol in dry-ice) and thawing (in water at room temperature with shaking) should be done rapidly.

Data Link UniProt <u>P03367</u> (gag-pol), UniProt <u>Q9YQ30</u> (HIV-1 Protease)

References : The HIV-1 strain and immunogen has been described in the following references.

- Adachi A *et al* "Production of acquired immunodeficiency syndrome-associated retrovirus in human and nonhuman cells transfected with an infectious molecular clone" J Virol 59: 284 -291(1986) PMID: <u>3016298</u>
- 2. Saitoh A *et a*l "Overproduction of human immunodeficiency virus type I reverse transcriptase in Escherichia coli and purification of the enzyme" *Microbiol Immunol* 34:509-521 (1990) PMID: <u>1699113</u>






Fig.2. Dot blotting of HIV-1 protease by using anti-HIV-1 protease antibody.

Anti-HIV-1 protease antibody was used at 1/2,000 dilution. As second antibody, goat anti-rabbit IgG antibody conjugated with HRP was used at 1/5,000 dilution. ECL system was used.

