

Cherwell Innovation Centre Unit 77, Heyford Park Upper Heyford Oxfordshire , OX25 5HD United Kingdom Tel: +44 (0)1869 238 067

Fax: +44 (0)1869 238 001

Certificate of Analysis

Name: Clostridium difficile Toxin A from ribotype 078

Product code: R78-TNA-10

Batch #: 9103

Description: Native toxin A purified from C. difficile strain 210708 ribotype 078

Purity: >90% pure by SDS-PAGE

Presentation: Lyophilised

Resuspension instructions

Contents when reconstituted in 109 μ l sterile distilled water will contain 10.028 μ g toxin A in 0.05M Hepes pH 7.4, 0.15M NaCl and 5% sucrose at a concentration of 92 μ g/ml

Usage guidelines

Storage: Lyophilised vial: +2°C to +8°C

Reconstituted liquid: +2°C to +8°C

Stability: Lyophilised vial: 1 year after date of receipt

Reconstituted liquid: 1 month

Freezing: Reconstituted liquid can be frozen at -80°C to extend stability to 1 year; however every

freeze/thaw cycle will cause increased aggregation. Avoid multiple freeze/thaw cycles and

multiple temperature fluctuations.

Mixing: Swirl or pipette gently to mix. Vortexing promotes aggregation.

Societ peruz & 11-500-2011

www.thenativeantigencompany.com VAT No: 102038475

contact@thenativeantigencompany.com Registration No. 7386339

