

## Data Sheet

### FluoTag-X2 anti-Chicken IgY Fluorophore-conjugated

Fluorophore	Catalog No. (20µl)	Catalog No. (200µl)
Atto488	N0702-At488-S	N0702-At488-L
Atto542	N0702-At542-S	N0702-At542-L
Atto565	N0702-At565-S	N0702-At565-L
Atto647N	N0702-At647N-S	N0702-At647N-L
Abberior-Star 580	N0702-Ab580-S	N0702-Ab580-L
Abberior-Star 635P	N0702-Ab635P-S	N0702-Ab635P-L
Sulfo-Cyanine 3	N0702-SC3-S	N0702-SC3-L
Sulfo-Cyanine 5	N0702-SC5-S	N0702-SC5-L

**Identity** Camelid sdAb anti-Chicken IgY

**Content** 5 µM fluorescently labeled sdAb in buffered saline, 50% glycerol, 0.09% sodium azide.

**Storage** Up to 3 months: -20 °C  
Up to 12 months: -80 °C or below  
Protect from light!

**Applications** IF / ICC      Recommended dilution 1:500

**Specificity** Recognizes specifically Chicken IgY  
Does not cross-react with antibodies from species like: mouse, rabbit, guinea pig, rat or goat.

For more information please visit our web page at [www.nano-tag.com](http://www.nano-tag.com)

To be used *in vitro* / for research only.  
Non-toxic, non-hazardous, non-infectious.