

MAP 2

Cat.No. 188 006; Polyclonal chicken antibody, 50 µg specific antibody (lyophilized)

Data Sheet

Reconstitution/ Storage	50 µg purified IgY, lyophilized. Ovalbumin was added for stabilization. For reconstitution add 50 µl H ₂ O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C until use. Before storing at -20°C add 1 vol of glycerol.
Applications	WB: 1 : 1000 (AP staining) IP: not tested yet ICC: 1 : 1000 IHC: 1 : 500 IHC-P/FFPE: 1 : 500
Immunogen	Recombinant protein corresponding to AA 2 to 314 from human MAP2-4 hu (UniProt Id: P11137-4)
Reactivity	Reacts with: rat (P15146), mouse (P20357). Other species not tested yet.
Specificity	Specific for MAP 2; recognizes all four isoforms.
matching control	188-OP

TO BE USED IN VITRO / FOR RESEARCH ONLY
NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

There are two major classes of heat stable microtubule associated proteins (MAPs): **MAP 2**, and tau. Both protein classes are involved in the regulation of microtubule polymerization in cells. Four differentially regulated isoforms of MAP 2 have been described so far.

Selected References SYSY Antibodies

Improved Methods for Fluorescence Microscopy Detection of Macromolecules at the Axon Initial Segment.
Alshammari MA, Alshammari TK, Laezza F
Frontiers in cellular neuroscience (2016) 10: 5. IHC

Selected General References

The MAP2/Tau family of microtubule-associated proteins.
Dehmelt L, Halpain S
Genome biology (2005) 6(1): 204.

MAP2 and tau bind longitudinally along the outer ridges of microtubule protofilaments.
Al-Bassam J, Ozer RS, Safer D, Halpain S, Milligan RA
The Journal of cell biology (2002) 157(7): 1187-96.

Distribution of map2 in hippocampus and cerebellum of young and old rats by quantitative immunohistochemistry.
Di Stefano G, Casoli T, Fattoretti P, Gracciotti N, Solazzi M, Bertoni-Freddari C
The journal of histochemistry and cytochemistry : official journal of the Histochemistry Society (2001) 49(8): 1065-6.

Age-dependent organotypic expression of microtubule-associated proteins (MAP1, MAP2, and MAP5) in rat brain.
Chauhan N, Siegel G
Neurochemical research (1997) 22(6): 713-9.

Variations in in vivo phosphorylation at the proline-rich domain of the microtubule-associated protein 2 (MAP2) during rat brain development.
Sánchez C, Díaz-Nido J, Avila J
The Biochemical journal (1995) 306 (Pt 2): 481-7.

Microtubule-associated proteins: a monoclonal antibody to MAP2 binds to differentiated neurons.
Izant JG, McIntosh JR
Proceedings of the National Academy of Sciences of the United States of America (1980) 77(8): 4741-5.