

Parvalbumin

Cat.No. 195 004; Polyclonal Guinea pig antibody, 100 µl antiserum (lyophilized)

Data Sheet

| | |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reconstitution/ Storage | 100 µl antiserum, lyophilized. For reconstitution add 100 µl H ₂ O, then aliquot and store at -20°C until use. |
| Applications | WB: 1 : 1000 (AP staining) IP: yes ICC: not tested yet IHC: 1 : 200 up to 1 : 500 IHC-P/FFPE: 1 : 500 EM: yes |
| Immunogen | Recombinant protein corresponding to AA 1 to 133 from rat Parvalbumin |
| Reactivity | Reacts with: rat (P02625), mouse (P32848), chicken, zebrafish, sheep, human (P20472). Other species not tested yet. |
| Specificity | Specific for parvalbumin. |
| matching control | 195-OP |
| Remarks | Due to its small size a tricine gel is recommended. |

TO BE USED IN VITRO / FOR RESEARCH ONLY NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

Parvalbumin is a small, acidic, calcium binding protein and belongs to the family of EF hand proteins. The protein is found in skeletal muscle and the brain of vertebrates where it locates to a specific population of GABAergic interneurons. This subset of neurons may contribute to maintaining the balance between excitation and inhibition in the cortex and the hippocampus.

Rudolf-Wissell-Str. 28
37079 Göttingen, Germany
Phone: +49 551-50556-0
Fax: +49 551-50556-384
E-mail: sales@sysy.com
Web: www.sysy.com

Selected References SYSY Antibodies

- A Hippocampus-Accumbens Tripartite Neuronal Motif Guides Appetitive Memory in Space.
Trouche S, Koren V, Doig NM, Ellender TJ, El-Gaby M, Lopes-Dos-Santos V, Reeve HM, Perestenko PV, Garas FN, Magill PJ, Sharott A, et al.
Cell (2019) : . **IHC, EM; tested species: mouse**
- Downregulation of Npas4 in parvalbumin interneurons and cognitive deficits after neonatal NMDA receptor blockade: relevance for schizophrenia.
Shepard R, Heslin K, Hagerdorn P, Coutellier L
Translational psychiatry (2019) 9(1): 99. **IHC; tested species: mouse**
- The synaptic ribbon is critical for sound encoding at high rates and with temporal precision.
Jean P, Lopez de la Morena D, Michalski S, Jaime Tobón LM, Chakrabarti R, Picher MM, Neef J, Jung S, Gültas M, Maxeiner S, Neef A, et al.
eLife (2018) 7: . **IHC; tested species: mouse**
- Neurochemical Characterization of PSA-NCAM+ Cells in the Human Brain and Phenotypic Quantification in Alzheimer's Disease Entorhinal Cortex.
Murray HC, Swanson MEV, Dieriks BV, Turner C, Faull RLM, Curtiss MA
Neuroscience (2018) 372: 289-303. **IHC; tested species: human**
- Supramammillary Nucleus Afferents to the Dentate Gyrus Co-release Glutamate and GABA and Potentiate Granule Cell Output.
Hashimoto-dani Y, Karube F, Yanagawa Y, Fujiyama F, Kano M
Cell reports (2018) 25(10): 2704-2715.e4. **IHC; tested species: mouse**
- Somatostatin Interneurons Facilitate Hippocampal-Prefrontal Synchrony and Prefrontal Spatial Encoding.
Abbas AI, Sundiang MJM, Henoch B, Morton MP, Bolkan SS, Park AJ, Harris AZ, Kellendonk C, Gordon JA
Neuron (2018) : . **IHC; tested species: mouse**
- Hippocampal CB1 Receptors Control Incidental Associations.
Busquets-Garcia A, Oliveira da Cruz JF, Terral G, Zottola ACP, Soria-Gómez E, Contini A, Martin H, Redon B, Varilh M, Ioannidou C, Drago F, et al.
Neuron (2018) : . **IHC; tested species: mouse**
- Automated analysis of images for molecular quantification in immunohistochemistry.
Guirado R, Carceller H, Castillo-Gómez E, Castrén E, Nacher J
Heliyon (2018) 4(6): e00669. **IHC; tested species: mouse**
- VIP-immunoreactive interneurons within circuits of the mouse basolateral amygdala.
Rhomberg T, Rovira-Esteban L, Vikár A, Paradiso E, Kremsner C, Nagy-Pál P, Papp OI, Tasan R, Erdélyi F, Szabó G, Ferraguti F, et al.
The Journal of neuroscience : the official journal of the Society for Neuroscience (2018) : . **IHC; tested species: mouse**
- High frequency neural spiking and auditory signaling by ultrafast red-shifted optogenetics.
Mager T, Lopez de la Morena D, Senn V, Schlotte J, D'Errico A, Feldbauer K, Wrobel C, Jung S, Bodensiek K, Rankovic V, Browne L, et al.
Nature communications (2018) 9(1): 1750. **IHC; tested species: mouse**
- Non-canonical heterogeneous cellular distribution and co-localization of CaMKIIα and CaMKIIβ in the spinal superficial dorsal horn.
Larsson M
Brain structure & function (2018) 223(3): 1437-1457. **IHC; tested species: rat**
- Shared rhythmic subcortical GABAergic input to the entorhinal cortex and presubiculum.
Viney TJ, Salib M, Joshi A, Unal G, Berry N, Somogyi P
eLife (2018) 7: . **IHC; tested species: mouse**
- Phase-specific surround suppression in mouse primary visual cortex correlates with figure detection behavior based on phase discontinuity.
Li F, Jiang W, Wang TY, Xie T, Yao H
Neuroscience (2018) : . **IHC; tested species: mouse**
- GABA and GABABreceptor subunit localization on neurochemically identified neurons of the human subthalamic nucleus.
Wu XH, Song JJ, Faull RLM, Waldvogel HJ
The Journal of comparative neurology (2018) 526(5): 803-823. **IHC; tested species: human**
- Parvalbumin-producing striatal interneurons receive excitatory inputs onto proximal dendrites from the motor thalamus in male mice.
Nakano Y, Karube F, Hirai Y, Kobayashi K, Hioki H, Okamoto S, Kameda H, Fujiyama F
Journal of neuroscience research (2018) : . **IHC; tested species: mouse**
- Ca²⁺-binding protein 2 inhibits Ca²⁺-channel inactivation in mouse inner hair cells.
Picher MM, Gehrt A, Meese S, Ivanovic A, Predeoehl F, Jung S, Schrauwen I, Dragonetti AG, Colombo R, Van Camp G, Strenzke N, et al.
Proceedings of the National Academy of Sciences of the United States of America (2017) 114(9): E1717-E1726. **IHC; tested species: mouse**