

SNAP 23

Cat.No. 111 203; Polyclonal rabbit antibody, 50 µg specific antibody (lyophilized)

Data Sheet

Reconstitution/ Storage	50 µg specific antibody, lyophilized. Affinity purified with the immunogen. Rabbit serum albumin was added for stabilization. For reconstitution add 50 µl H ₂ O to get a 1mg/ml solution in PBS. Then aliquot and store at -20°C until use.
Applications	WB: 1 : 1000 (AP staining) IP: yes (see remarks) ICC: 1 : 50 IHC: yes IHC-P/FFPE: not tested yet
Immunogen	Synthetic peptide corresponding to AA 196 to 211 from human SNAP23 (UniProt Id: O00161)
Reactivity	Reacts with: human (O00161), rat (O70377), mouse (O09044), hamster, pig, rabbit. Other species not tested yet.
Specificity	Specific for SNAP 23. (K.D. verified)
matching control	111-2P
Remarks	IP: Antibody 2 (cat. no. 111 213) is recommended.

TO BE USED IN VITRO / FOR RESEARCH ONLY NOT TOXIC, NOT HAZARDOUS, NOT INFECTIOUS, NOT CONTAGIOUS

SNAP 23 (synaptosome-associated protein of 23 kDa) is a ubiquitously expressed isoform and functional homologue of SNAP 25. It is resistant to cleavage by BoNT/A and E. The protein is part of the exocytotic fusion complex (v-SNARE) where it assembles with syntaxin 1 and synaptobrevin. SNAP 23 is able to function in regulated exocytosis. Both isoforms may have their own specific binding partners and discrete, albeit mechanistically similar, functional roles within the cell.

Selected References SYSY Antibodies

- Cholesterol-dependent syntaxin-4 and SNAP-23 clustering regulates caveolar fusion with the endothelial plasma membrane. Predescu SA, Predescu DN, Shimizu K, Klein IK, Malik AB
The Journal of biological chemistry (2005) 280(44): 37130-8. **WB, ICC**
- Identification of soluble N-ethylmaleimide-sensitive factor attachment protein receptor exocytotic machinery in human plasma cells: SNAP-23 is essential for antibody secretion. Reales E, Mora-López F, Rivas V, García-Poley A, Brieva JA, Campos-Caro A
Journal of immunology (Baltimore, Md. : 1950) (2005) 175(10): 6686-93. **WB, ICC**
- Distribution of synaptic vesicle proteins in the mammalian retina identifies obligatory and facultative components of ribbon synapses. Von Kriegstein K, Schmitz F, Link E, Südhof TC
The European journal of neuroscience (1999) 11(4): 1335-48. **WB, ICC**
- SNARE-dependent glutamate release in megakaryocytes. Thompson CJ, Schilling T, Howard MR, Genever PG
Experimental hematology (2010) 38(6): 504-15. **IHC**
- Phosphorylation of syntaxin-3 at Thr 14 negatively regulates exocytosis in RBL-2H3 mast cells. Tadokoro S, Shibata T, Inoh Y, Amano T, Nakanishi M, Hirashima N, Utsunomiya-Tate N
Cell biology international (2016) 40(5): 589-96. **WB; tested species: rat**
- Syntaxin 8 regulates platelet dense granule secretion, aggregation, and thrombus stability. Golebiewska EM, Harper MT, Williams CM, Savage JS, Goggs R, Fischer von Mollard G, Poole AW
The Journal of biological chemistry (2015) 290(3): 1536-45. **WB**
- Insulin stimulates syntaxin4 SNARE complex assembly via a novel regulatory mechanism. Kioumourtzoglou D, Gould GW, Bryant NJ
Molecular and cellular biology (2014) 34(7): 1271-9. **ICC**
- Renal epithelial cells can release ATP by vesicular fusion. Bjaelde RG, Arnadottir SS, Overgaard MT, Leipziger J, Praetorius HA
Frontiers in physiology (2013) 4: 238. **WB; KD verified; tested species: rabbit**
- REST/NRSF governs the expression of dense-core vesicle gliosecretion in astrocytes. Prada I, Marchaland J, Podini P, Magrassi L, D'Alessandro R, Bezzi P, Meldolesi J
The Journal of cell biology (2011) 193(3): 537-49. **WB**
- A role for VAMP8/endobrevin in surface deployment of the water channel aquaporin 2. Wang CC, Ng CP, Shi H, Liew HC, Guo K, Zeng Q, Hong W
Molecular and cellular biology (2010) 30(1): 333-43. **WB**
- Promiscuous interaction of SNAP-25 with all plasma membrane syntaxins in a neuroendocrine cell. Bajohrs M, Darios F, Peak-Chew SY, Davletov B
The Biochemical journal (2005) 392(Pt 2): 283-9. **ICC**
- Syntaxin 1A has a specific binding site in the H3 domain that is critical for targeting of H⁺-ATPase to apical membrane of renal epithelial cells. Li G, Yang Q, Alexander EA, Schwartz JH
American journal of physiology. Cell physiology (2005) 289(3): C665-72. **WB**
- PKA-dependent ENaC trafficking requires the SNARE-binding protein complexin. Butterworth MB, Frizzell RA, Johnson JP, Peters KW, Edinger RS
American journal of physiology. Renal physiology (2005) 289(5): F969-77. **ICC**
- Vesicular glutamate transporter-dependent glutamate release from astrocytes. Montana V, Ni Y, Sunjara V, Hua X, Parpura V
The Journal of neuroscience : the official journal of the Society for Neuroscience (2004) 24(11): 2633-42. **WB**
- Intersectin regulates fission and internalization of caveolae in endothelial cells. Predescu SA, Predescu DN, Timblin BK, Stan RV, Malik AB
Molecular biology of the cell (2003) 14(12): 4997-5010. **WB**
- Syntaxin isoform specificity in the regulation of renal H⁺-ATPase exocytosis. Li G, Alexander EA, Schwartz JH
The Journal of biological chemistry (2003) 278(22): 19791-7. **WB**
- Vimentin filaments in fibroblasts are a reservoir for SNAP23, a component of the membrane fusion machinery. Faigle W, Colucci-Guyon E, Louvard D, Amigorena S, Galli T
Molecular biology of the cell (2000) 11(10): 3485-94. **ICC**
- EHS1/intersectin, a protein that contains EH and SH3 domains and binds to dynamin and SNAP-25. A protein connection between exocytosis and endocytosis? Okamoto M, Schoch S, Südhof TC
The Journal of biological chemistry (1999) 274(26): 18446-54. **WB**